





OPTIMA CL323i

NOTE: THIS IS A SUGGESTED LAYOUT BASED OFF OF GE HEALTHCARE PRE-INSTALL MANUAL. PLEASE CONTACT LOCAL PROJECT MANAGER OF INSTALLATION IF YOUR ROOM LAYOUT OR DIMENSIONS ARE DIFFERENT THAN THE LAYOUT SHOWN.



NOTE: THIS DRAWING TO BE SCALED UP BY TWO.

ROOM DIMENSIONS	LENGTH & WIDTH	CEILING HEIGHT
RECOMMENDED:	32'-0" x 20'-2" [9.7m x 6.1m]	9'-6" [2.9m]



SEE APPENDIX "A" FOR POWER SPECIFICATIONS

ROOM PLAN LOCATED ON PAGE 4-79

ITEM	QUANTITY	ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT
①	1	ATLAS CABINET <C2>	659 lbs	1825 btu
②	1	ATLAS CABINET <C1>	1115 lbs	3389 btu
③	1	COUNTERBALANCED EYE AND THYROID SHIELD WITH R96 LAMP	143 lbs	
④	1	AW WORKSTATION	81 lbs	1201 btu
⑤	1	FOUR LCD MONITOR SUSPENSION ON 9 ft. 6 in. XT INBOARD BRIDGE	485 lbs	1157 btu
⑥	2	LONGITUDINAL STATIONARY RAIL FOR XT SUSPENSION	68 lbs	
⑦	1	INNOVA POSITIONER <REFERENCE TABLE BASE-PLATE DETAIL FOR FLOOR MOUNTING INFORMATION>	1653 lbs	2416 btu
⑧	1	TRAM NET RACK	8 lbs	
⑨	1	COLOR PRINTER		1054 btu
⑩	1	REMOTE CONTROL FOR INJECTOR	4 lbs	
⑪	1	INJECTOR HEAD ON TABLE RAIL	15 lbs	
⑫	1	INJECTOR ELECTRONICS ON SHELF	37 lbs	320 btu
⑬	1	OMEGA V TABLE	1300 lbs	614 btu
⑭	1	XR BUZZER <LOCATED ABOVE CEILING>	2 lbs	
⑮	1	UPS INTERFACE BOX		
⑯	1	UPS CABINET	1170 lbs	4061 btu
⑰	1	DETECTOR CHILLER	33 lbs	706 btu
⑱	1	WATER CHILLER	449 lbs	18716 btu
⑲	1	CONTROL ROOM MONITOR WITH DL KEYPAD	22 lbs	204 btu
⑳	1	OPERATORS CONSOLE	22 lbs	546 btu
㉑	1	BOLUS CHASE HANDSWITCH	2 lbs	
㉒	1	WORKSTATION CART		
㉓	1	MAC-LAB CONSOLE, INCLUDES MONITORS AND KEYBOARD	566 lbs	2935 btu
⑤①	1	VITALINQ MICROPHONE <ONE ON MONITOR BRIDGE IN EXAM ROOM>		
⑤②	1	VITALINQ MICROPHONE <ONE ON COUNTERTOP IN CONTROL ROOM>		
⑤③	6	VITALINQ SPEAKER		
⑤④	1	INNOVA MAIN DISCONNECT, REFERENCE JUNCTION POINT 'PDB' ON SHEET E1 FOR DETAILED DESCRIPTION.	326 lbs	1532 btu
⑤⑤	1	VITALINQ CONSOLE		

ENVIRONMENT

Ambient Operation Temperature: 68°-77° F, (20°-25° C)
 Allowable Temperature Change: Not Available
 Humidity: 30%-75% non condensing
 Allowable Humidity Change: Not Available
 Max - Altitude: 9842 ft. (3000 M) above sea level

ROOM PLAN LOCATED ON PAGE 4-79

ITEM	QUANTITY	ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT
		***** CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS *****		
60		COUNTER TOP FOR EQUIPMENT-SHELVING MAY BE REQUIRED PROVIDE GROMMETED OPENINGS AS REQUIRED TO ROUTE INTERCONNECT CABLES TO RACEWAY BELOW COUNTERTOP.		
61		CONTROL WALL TO CEIL. WITH LEAD GLASS VIEWING WINDOW.		
62		COUNTER TOP WITH SINK, BASE AND WALL CABINETS.		
63		SCRUB SINK		
64		COUNTER TOP WITH BASE AND WALL CABINETS		
65		BEARING BLOCK OUTLINE, SEE S1 FOR MORE INFORMATION.		
66		CABLE DRAPE RAIL.		
67		MED GASES IN CEILING		
68		CUSTOMER SUPPLIED		
69		CATHETER CABINETS		
70		LEAD APRON RACK		
71		SHELF - CUSTOMER TO PROVIDE ADEQUATE WALL SUPPORT		
72		X-RAY ON WARNING LIGHT - AVAILABLE FROM GE SUPPLY CALL 800-200-9760 GE CAT. NO. WXIABWW-DF-XIU		
73		MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 44 IN. W x 83 IN. H [1118mm x 2108mm], CONTINGENT ON A 96 IN. [2438mm] CORRIDOR WIDTH		
74		150-AMP LOCAL SERVICE DISCONNECT FOR LOCK-OUT/ TAG-OUT CAPABILITY. (MAY BE A FUSED DISCONNECT, CIRCUIT BREAKER OR SAFETY SWITCH.)		
90		X-RAY ROOM WARNING LIGHT/ROOM LIGHTING CONTROL PANEL REFERENCE JUNCTION POINT 'XRLC' ON SHEET 'E1' FOR DETAILED DESCRIPTION -CAT. NO. E4502SS FOR WARNING LIGHT & ROOM LIGHT CONTROL.		