





INNOVA IGS 520, 530, 540 INTERVENTIONAL RADIOLOGY (IR)



NOTE: THIS IS A SUGGESTED LAYOUT BASED OFF OF GE HEALTHCARE PRE-INSTALL MANUAL. PLEASE CONTACT LOCAL PROJECT MANAGER OF INSTALLATION IF YOUR ROOM LAYOUT OR DIMENSIONS ARE DIFFERENT THAN THE LAYOUT SHOWN.

ROOM DIMENSIONS	LENGTH & WIDTH	CEILING HEIGHT
RECOMMENDED:	32'-0" x 25'-0" [9.75m x 7.62m]	9'-6" [2.89m]



SEE APPENDIX "A" FOR POWER SPECIFICATIONS

ROOM PLAN LOCATED ON PAGE 4-67

ITEM	QUANTITY	ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT
①	1	INJECTOR ELECTRONICS	37 lbs	320 btu
②	1	TRAM NET RACK	8 lbs	
③	1	XR BUZZER (LOCATED ABOVE CEILING)	2 lbs	
④	1	INNOVA IQ TABLE	1750 lbs	614 btu
⑤	1	INNOVA POSITIONER (REFERENCE TABLE BASE-PLATE DETAIL FOR FLOOR MOUNTING INFORMATION)	1653 lbs	2416 btu
⑥	1	UPS INTERFACE BOX		
⑦	1	ATLAS CABINET (C2)	659 lbs	1825 btu
⑧	1	ATLAS CABINET (C1)	1115 lbs	3389 btu
⑨	1	UPS CABINET	1170 lbs	4061 btu
⑩	1	DETECTOR CHILLER	33 lbs	706 btu
⑪	1	COOLIX 4100 AUTOTRANSFORMER	66 lbs	238 btu
⑫	1	COOLIX 4100 WATER CHILLER	264 lbs	18730 btu
⑬	1	IVUS VOLCANO S51 CONSOLE, INCLUDES FLAT PANEL MONITOR AND KEYBOARD (DESK MOUNTED)	68 lbs	1631 btu
⑭	1	LARGE DISPLAY MONITOR ON SINGLE MONITOR SUSPENSION 9 ft. 6 in. INBOARD BRIDGE (MOUNT TWO GE MONITORS ON BACKSIDE OF LD MONITOR)	784 lbs	1706 btu
⑮	1	LARGE DISPLAY MONITOR CABINET	253 lbs	3412 btu
⑯	1	3 KVA UPS CABINET (LARGE DISPLAY SUBSYSTEM OPTION)	77 lbs	546 btu
⑰	1	19 in. MONITOR ON WALL SUPPORT	26 lbs	204 btu
⑱	1	CONTROL ROOM MONITOR WITH DL KEYPAD	22 lbs	204 btu
⑲	1	OPERATORS CONSOLE	22 lbs	546 btu
⑳	1	AW WORKSTATION	81 lbs	1201 btu
㉑	1	VIVID 1 ULTRASOUND UNIT (ON SAFELOCK CART)		
㉒	1	NURSING NOTES WORKSTATION	46 lbs	682 btu
㉓	1	COUNTERBALANCED EYE AND THYROID SHIELD WITH R96 LAMP	143 lbs	
㉔	2	LONGITUDINAL STATIONARY RAIL FOR XT SUSPENSION	68 lbs	
㉕	1	REMOTE CONTROL FOR INJECTOR	4 lbs	
㉖	1	INJECTOR HEAD ON TABLE RAIL	15 lbs	

ENVIRONMENT

Ambient Operation Temperature: 68°-77° F, (20°-25° C)
 Allowable Temperature Change: Not Available
 Humidity: 30%-75% non condensing
 Allowable Humidity Change: Not Available
 Max - Altitude: 9842 ft. (3000 M) above sea level

ROOM PLAN LOCATED ON PAGE 4-67

ITEM ⬡	QUANTITY	ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT
27	1	IVUS VOLCANO COLOR PRINTER		
28	1	BOLUS CHASE HANDSWITCH	2 lbs	
29	1	WORKSTATION CART		
30	1	MAC-LAB CONSOLE, INCLUDES MONITORS AND KEYBOARD	566 lbs	2935 btu
50	1	INNOVA MAIN DISCONNECT, REFERENCE JUNCTION POINT 'PDB' ON SHEET E1 FOR DETAILED DESCRIPTION.	326 lbs	1532 btu
51	1	VITALINQ CONSOLE		
52	6	VITALINQ SPEAKER		
53	1	VITALINQ MICROPHONE		
54	1	VITALINQ MICROPHONE (ONE ON MONITOR BRIDGE IN EXAM ROOM)		
***** CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS *****				
60		X-RAY ON WARNING LIGHT - AVAILABLE FROM GE SUPPLY CALL: 800-200-9760 GE CAT. NO. WXIABWW-OF-XIU		
61		MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 44 IN. W x 83 IN. H [1118mm x 2108mm], CONTINGENT ON A 96 IN. [2438mm] CORRIDOR WIDTH		
62		150-AMP LOCAL SERVICE DISCONNECT FOR LOCK-OUT/TAG-OUT CAPABILITY. (MAY BE A FUSED DISCONNECT, CIRCUIT BREAKER OR SAFETY SWITCH.)		
63		CATHETER CABINETS		
64		CUSTOMER SUPPLIED STORAGE CABINET		
65		SHELF - CUSTOMER TO PROVIDE ADEQUATE WALL SUPPORT		
66		BEARING BLOCK OUTLINE, SEE S1 FOR MORE INFORMATION.		
67		CABLE DRAPE RAIL.		
68		MED GASES IN CEILING		
69		LEAD APRON RACK		
70		COUNTER TOP FOR EQUIPMENT- MINIMUM DEPTH 30 in. OR ADDITIONAL SHELVING MAY BE REQUIRED PROVIDE GROMMETED OPENINGS AS REQUIRED TO ROUTE INTERCONNECT CABLES TO RACEWAY BELOW COUNTERTOP.		
71		FLUSH ACCESS FLOORING		
90		X-RAY ROOM WARNING LIGHT/ROOM LIGHTING CONTROL PANEL REFERENCE JUNCTION POINT 'XRLC' ON SHEET 'E1' FOR DETAILED DESCRIPTION -CAT. NO. E4502SS FOR WARNING LIGHT & ROOM LIGHT CONTROL.		