





Lunar iDXA System



ROOM DIMENSIONS	LENGTH & WIDTH	CEILING HEIGHT
MINIMUM:	11'-0" x 13'-2" [3.35m x 4.0m]	7'-0" [2.1m]



SEE APPENDIX "A" FOR POWER SPECIFICATIONS

ROOM PLAN LOCATED ON PAGE B-2

ITEM	QUANTITY	ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT
①	1	WORKSTATION	28 lbs	399 btu
②	1	LUNAR IDXA SCAN TABLE	793 lbs	150 btu
<p>***** CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS *****</p>				
90		1 KVA ISOLATION TRANSFORMER; GE PART LU8544		
91		SINGLE PHASE CIRCUIT BREAKER WITH ENCLOSURE, GE ENERGY PART #TLM612FCUD, THQL1115 & TGK12 AVAILABLE FROM GEXPRO MILWAUKEE 414-527-6600		
92		OPERATORS CHAIR		

ENVIRONMENT

Ambient Operation Temperature: 65-81° F (18-27° C)
 Allowable Temperature Change: 3.6° F (1.5° C) per hour
 Humidity: 20-80 percent non condensing
 Allowable Humidity Change: 10 percent per hour
 Max - Altitude: Undefined in Specifications