

ROOM PLAN LOCATED ON PAGE 8-197

ITEM	QUANTITY	ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT
①	1	SYSTEM CABINET	2698 lbs	9556 btu
②	1	MESH SHIELD		
③	1	PENETRATION PANEL		
④	1	MAGNET RUNDOWN UNIT	8 lbs	
⑤	1	OPERATOR WORKSPACE W/COLOR LCD MONITOR	125 lbs	4948 btu
⑥	1	OPERATOR'S CHAIR		
⑦	1	OPERATOR WORKSPACE CABINET	198 lbs	
⑧	1	PATIENT ALERT CONTROL BOX		
⑨	1	1.5 TESLA LCC ACTIVE SHIELD MAGNET	13115 lbs	8191 btu
⑩	1	PATIENT TRANSPORT TABLE (DOES NOT INCLUDE PATIENT)	279 lbs	
⑪	1	COOLING UNIT	121 lbs	5699 btu
⑫	1	BLOWER BOX	46 lbs	3412 btu
⑬	1	ABSORBER	99 lbs	
⑭	1	MAGNET MONITOR	22 lbs	204 btu
⑮	1	SPT PHANTOM CABINET	350 lbs	
⑯	0	INJECTOR HEAD ON PEDESTAL (OPTIONAL)	59 lbs	
⑰	0	CONTROL ROOM UNIT (OPTIONAL)	15 lbs	
⑱	0	BATTERY CHARGING UNIT (OPTIONAL)	4 lbs	
⑵	1	MR COMMON CHILLER SYSTEM	683 lbs	
⑶	1	SHIELD COOLER COMPRESSOR	253 lbs	31379 btu
⑷	1	REMOTE CONTROL FOR CHILLER SYSTEM	2 lbs	

ENVIRONMENT

Ambient Operation Temperature: 50-70° F [10-21° C]
 Allowable Temperature Change: 5° F [3° C] per hour
 Humidity: 30-60 per cent non condensing
 Allowable Humidity Change: 10 per cent per hour
 Max - Altitude: 8000 ft. [2438m] above sea level

ROOM PLAN LOCATED ON PAGE 8-197

ITEM	QUANTITY	ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT
		***** CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS *****		
[60]		RF FILTERS - LOCATE WITHIN 24 in. RF COMMON GROUND STUD [610 mm] OF THE		
[61]		AIR CONDITIONING. <VIBRATION ISOLATION IS RECOMMENDED AT SUPPORTS OF EACH UNIT TO BE INSTALLED.>		
[62]		MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 43 IN. W x 82 IN. H [1092mm x 2083mm], CONTINGENT ON A 96 IN. [2438mm] CORRIDOR WIDTH		
[63]		SHELF - CUSTOMER TO PROVIDE ADEQUATE WALL SUPPORT		
[64]		VALVES AND HOSE BARBS FOR COOLING SYSTEM		
[65]		MAGNET ROOM EXHAUST FAN		
[66]		MINIMUM 9 FT. -0 IN. [2743 mm] x 9 FT. -0 IN. [2743 mm] REMOVABLE WALL SECTION FOR MAGNET DELIVERY/REMOVAL.		
[67]		RF SCREEN, INCLUSIVE OF WALLS, FLOOR, DOOR, ETC. GROUND IMPEDANCE GREATER THAN 1000 OHMS. ATTENUATION 100dB AT 10-100MHZ PLANEWAVE.		
[68]		BASE CABINET FOR STORAGE OF: SURFACE COILS, PATIENT POSITIONING PADS, PHANTOMS, ETC.		
[90]		DC LIGHTING CONTROL PANEL 155 lbs <70 kg> 1024 BTU/HR. <CAT. NO. E4502SC/SE - BASIC SYSTEM>		
[91]		MAIN DISCONNECT CONTROL CAT NO. M3088TM		
[92]		DC LIGHTING AUTO TRANSFORMER 60 lbs [27 kg] <PART OF VARIABLE DIMMER SYSTEM> <CAT. NO. E4502SD/SF INCLUDES BASIC SYSTEM>		
[93]		METAL DETECTOR <HAND HELD>		